Measurement Unit (MU3)

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CYBRES[®] Measurement Unit (MU3)


for electrochemical and electrophysiological analysis of fluids and organic tissues

- Differential Electrochemical Impedance Spectrometer (EIS)

- Phytosensing and phytoactuating system

- Biosensor based on fermentation activity of yeast

The MU3 is a bio-hybrid interface device for real-time interactions with different fluidic, biological and microbiological systems. It includes differential Electrochemical Impedance Spectrometer, analyzer of bio-potentials, different bio- and environmental sensors. MU3 is capable of performing high-resolution differential measurements, where ionic properties of two fluidic or organic samples are compared with each other. Based on this approach, the system can perform (electro-)physiological analysis of plants and microorganisms, measure dynamic or static properties of corresponding tissues and solutions. Such task appears in applications, where e.g. ultra-weak electrochemical or (electro-) physiological changes should be detected, caused, among others, by non-chemical methods, environment or different technologies. The system is developed for single-run measurements or for a long-term monitoring with online graphical output in web. The detectors-actuators (DA) module executes real-time data processing and decision making for operation of different actuators. This functionality is useful for performing fully autonomous experiments and development of complex feedback-based and adaptive scenarios with electrochemical, biological and bio-hybrid systems.

The structure of MU system is shown in Figure above. It consists of an active module with electronic components and replaceable passive electrodes for phytosensing, biosensing and EIS applications.

Main features

- main processor: ARM cortex M3 MPU, 80 MHz
- hardware support of analysis: PSoC system
- non-volatile (flash) memory: 512 Mb
- sampling frequency: (12-24 bits) up to 1 Msps
- frequency range of EIS: 8Hz-200kHz

- conductivity ranges: 0.6 μ S/cm-200 mS/cm
- ranges of excitation voltage AC: 0.001-0.01V, 0.01-0.1V, 0.1-1V
- amplification factors: 50, 500, 5000, 50000
- EIS analysis: Amplitudes, FRA Phase, RMS Magnitude, Correlation, Electrochemical stability in time/frequency/time-frequency domains, Statistical analysis, Excitation analysis
- measurement modes: 1) impedance spectrometer; 2) signal scope; 3) continuous measurements at a constant f ; 4) continuous measurements at variable f ; 5) Frequency Response Profile (FRP) at a fixed set of frequencies; 6) continuous FRP
- additional sensors: 3D accelerometer/magnetometer, two internal temp. sensors, external high-resolution temp. sensor, external high-resolution environmental data logger (optional)
- basic accuracy class: 0.5%, 0.1%
- powering: 5V (~0.3A), external active USB3.0 hub,
- weight/size: 112g, 96mm x 22mm x 70mm

Applications


1) General electrochemical applications are precise industrial fluidic measurements and differential fluidic meters in research and laboratory usage, detectors of weak (non-)electromagnetic emissions by analysing electrochemical changes in fluids. The device is suitable for the analysis of differential electrochemical properties of samples exposed by non-chemical, non-temperature, non-acoustic, non-mechanical and non-electromagnetic factors.

2) General electrophysiological applications are plant, organic tissue and microbiological measurements for phytosensing and biosensing usage, e.g. monitoring plant physiology and electrophysiology, analysis of bio-potentials and tissue conductivity, bio-sensors based on fermentation, sedimentation, gas production (or degassing), metabolic production or any other processes that change concentration and mobility of ions. The MU EIS system is designed for long-term monitoring of biological samples, e.g. for quality control purposes or for the analysis of biochemical reactions.

3) Analysis and measurements of weak interactions, in particular in research of certain quantum phenomena appearing in macroscopic systems. Examples are the proton tunneling effect and self-ionization of water based on quantum fluctuation of E-field. These quantum effects on the micro-level between water molecules, ions and protons, causes changes of fluidic parameters on the macro-level, which can be in turn measured as changes of e.g. impedance. The device allows statistically significant measurements of these effects with the standard EIS method.

Software installation (Windows 7, 8, 10)

1. Install the redistributable package for visual C++ 2012 (32/64 bit versions are in the directory 'drivers').
2. After connecting the device first time to PC, Windows will detect the MU EIS device as unrecognized 'USBUART' device. In the Device Manager with the right button click open the device settings and in manual mode show the pass to the 'USBUART_cdc.inf' file (file is included in the directory 'drivers'). After completing that the MU EIS device should appear as 'USB Serial Device (COM X)'.

3. The client program does not require installation. All the necessary files are contained in MU-EIS-Client directory.
4. Plotting is carried out by any program that can read numerical data from files (e.g. Microsoft Excel). The developers propose to use the free program gnuplot, however the decision to use this software lies entirely on users. To install the gnuplot program, run the installation file 'gp503-win64-mingw.exe' and install the program in the default directory. During the installation make sure that PATH environmental variable is set. If the gnuplot was installed into another directory, or the correct installation failed, you need to set the correct path and also in the PATH command. These steps need to perform only once when software is first time installed into your PC. Updates of the client program are performed only by replacing the program folder.

Getting started

1. Connect the USB cable to the device and select the COM port.
2. Configure the application (EIS, biosensor, phytosensor) and connect the selected electrodes.
3. For real-time plot, check the checkbox "online plot enable".
4. Press the button "Measurement Start".

Measurement modes and device versions

N	Application	DDS mode	Configur.
1	signal distortion analysis, fast statistical analysis	signal scope	EIS
2	analysis of temporal EIS/temperature dynamics, „experiment during measurement” mode, biosensor applications, 3D time-frequency analysis	continuous modes	EIS, biosensor
3	regression enabled , the highest resolution of temporal dynamics	continuous modes	EIS
4	excitation enabled , analysis of pre-treated fluidic samples, characterization of non-chemical treatment	continuous modes	EIS
5	frequency impedance spectroscopy, differential analysis of samples, FRA profiles	frequency modes	EIS
6	environmental measurements	off	EIS
7	electrophysiological measurements of tissues, physiological measurements (with corresponding sensors), differential potential analysis	off	phytosensor
8	electrochemical interface to bio-samples, biological tissues and plants	continuous modes	phytosensor, biosensor

N	Device versions	Hardware	Software
1	phytosensor basic	1 channel biopotentials & impedance	enabled
2	phytosensor advanced	+ TransAmb sensor (leaf transpiration) + 2 channels biopotentials & impedance	enabled
3	phytosensor full	+ sup flow sensor	enabled
4	EIS, open electrodes	+ 2x EIS electrodes, + fluidic/envir. t sensors, + excitation spectroscopy	activation code*
5	EIS, embedded version	different device with thermostat, + RGB/IR excitation spectroscopy	enabled
6	biosensor	+ fermentation module, + RGB/IR excitation spectroscopy	enabled

*transition from phytosensor to EIS/Biosensor requires software activation

All measurements in differential modes (time-differential or channel-differential) do not require calibration. The calibration can be performed by users for measuring an absolute value of conductivity at a fixed frequency or for linearizing the frequency dependent EIS dynamics. The device is CE/FCC/CCC/WEEE certified, certificates can be found on the device home page or by request: CYBRES GmbH, Melunerstr. 40, 70569 Stuttgart, Germany, info@cybertronica.de.com. Guarantee: 12 months after selling. Please register per email for obtaining software/hardware updates.

Documentation, user manual, application notes and other materials:

- www.cybertronica.de.com/products/MU-EIS-spectrometer
- www.cybertronica.de.com/products/phytosensor
- www.cybertronica.de.com/products/biosensor

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declare that this DoC is issued under our sole responsibility and belongs to the following product:

- product name: Measurement Unit (MU3)
- type: **CYBRES MU3**
- serial numbers: CYBRES MU3 34xxxx


The object of the declaration described above is in conformity with the relevant Community harmonisation legislation:

- The electromagnetic compatibility (EMC) Directive 2014/30/EU, test specifications: EN 61326-1:2013, IEC 61326-1:2012
- Restriction of the use of certain Hazardous Substances Directive (RoHS) 2011/65/EU, EN 50581 Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

Testing laboratory: SGS Germany GmbH, Hofmannstr. 50, D-81379 Munich, Test Report No.: N3P90001

Place, date: Stuttgart March 09, 2019

Signature:

Dr. Serge Kernbach, CEO

FCC Declaration of Conformity

The device **CYBRES MU3** complies with Part 15 of the FCC Rules and Regulations for Information Technology Equipment. Operation is subject to the following two conditions: (1) these devices may not cause harmful interference, and (2) these devices must accept any interference received, including interference that may cause undesired operation.

Product:

- product name: Measurement Unit (MU3)
- type: **CYBRES MU3**
- serial numbers: CYBRES MU3 34xxxx

Responsible Party:


- CYBRES GmbH, Melunerstr.40, 70569 Stuttgart, Germany
- Phone/Fax: +49-711-41001901
- info@cybertronica.de.com, www.cybertronica.de.com
- Contact person: Dr. Serge Kernbach

We, the responsible party CYBRES GmbH, declare that the products “Measurement Unit (**CYBRES MU3**)” was tested by the accredited laboratory SGS Germany GmbH to conform to the applicable FCC rules and regulations. The method of testing was in accordance with the appropriate measurement standards, and all necessary steps have been taken to ensure that all production units of these devices will continue to comply with the Federal Communications Commission's requirements.

- Test Specifications: [covered by accreditation] FCC 47 CFR Part 15, §15.107, §15.109, ICES -003 Issue 6.
- Accredited testing laboratory: SGS Germany GmbH, Hofmannstr. 50, D-81379 Munich, Test Report No.: N3P90002

Place, date: Stuttgart, March 09, 2019

Signature:

Dr. Serge Kernbach, CEO